

DRÄGER DRUGCHECK® 3000



The Dräger DrugCheck® 3000 Saliva Test Kit significantly reduces the issue of a false negative and positive results. The device has two indicators that inform the tester that they can proceed to the next step, confidently knowing that the procedure is correct.

The test kit rapidly tests individuals simultaneously for up to six substances, including Amphetamines, Benzodiazepines, Cannabis (THC), Cocaine, Methamphetamines and Opiates.



BENEFITS

- Find out within minutes if a person recently consumed certain drugs
- It does not require electricity and can be used anywhere
- The compact and quick oral fluid- based drug test yields reliable results affordably and efficiently
- High sensitivity (THC down to 15ng/ml)
Operating temperature 5 - 40°C



APPLICATIONS

- Workplace health and safety. mining, logistics, manufacturing, agriculture, among others
- Law enforcement
- Correctional facilities
- Hospitals and rehabilitation clinics
- Personal use



OUR SERVICES



Custom Systems - Data Loggers - Breathalysers
Drug Testing Kits - Post Harvest Technologies -
Environmental Monitoring

CONTACT US

+617 3361 2000 (Australia)



DRÄGER DRUGCHECK® 3000



Benefits & Features

The disposable **DrugCheck® 3000** cannot be manipulated and is hygienic to use.



Unambiguous on-the-spot drug detection

Its compact design makes the DrugCheck® 3000 easy to transport. It can be made ready for use quickly and easily and is suitable for roadside drug checks and drug testing in safety-relevant workplaces. The test kit has no electrical parts, making the drug tester safe to use in hazardous areas and even as point-of-care-testing (POCT).




Fast and sensitive analysis for up to six substance classes

You can check individuals for five substance classes simultaneously with Dräger DrugCheck® 3000: cocaine, opiates, amphetamine, methamphetamines/designer drugs (e.g. Ecstasy, MDMA), benzodiazepines (e.g. in prescription medicines) and cannabis (THC).



Dräger Drugcheck App and QR Scan

Before performing a test, the QR on the test can be scanned. User and subject details are entered, and test performed. Once results are visible, the app will scan the result, photograph it, and email it to your records' nominated location. Simple.



OUR SERVICES



Custom Systems - Data Loggers - Breathalysers
Drug Testing Kits - Post Harvest Technologies -
Environmental Monitoring

CONTACT US

+617 3361 2000 (Australia)



www.pacdatasys.com.au



sales@pacdatasys.com.au



27 High Tech Court, Eight
Mile Plains, QLD, 4113